

Side-by-Side



The Side-by-Side tricycle is much more than just a bike. It represents regained freedom, “wind in your hair”, rosy cheeks and vitality. It enriches the lives of people with walking difficulties or disabilities each day

The Side-by-Side design makes the tricycle ideal for the disabled and those with severe walking disabilities, either as an outdoor leisure activity or for structured rehabilitation.

The Side-by-Side is renowned for its user-friendliness, flexibility and excellent performance:

- easy to access and safe to ride
- comfortable and easy adjustable seats with removable seat cushions
- available with a powerful electric motor
- rear wheel drive via a differential on the rear axle
- extremely rugged with a long lifetime
- recognized rehabilitation aid
- extensive accessories list
- adaptable to your special needs and delivered within a few weeks

Standard equipment:

- 7 speed
- Differential
- Foot and hand brake
- Parking brake
- Adjustable seats and handlebars
- Front basket
- Lock with 3 keys
- Rear view mirror
- Seat cushions
- Schwalbe Marathon Plus tyres

Side-by-Side

Side-by-Side

The Side-by-Side is well suited for rehabilitation purposes for the elderly and for people with walking difficulties. It offers the function of linking the cadence of pilot and passenger together, so the driver stimulates the passenger to pedal along, which is naturally important in the rehabilitation process.

This feature can easily be turned off, offering the passenger the choice of pedaling or taking a rest, while the driver is still pedaling.



Side-by-side without motor

Side-by-Side P20

(250W front wheel motor)

The Side-by-Side P20 is equipped with a speed sensor that measures pedal use, which is key for getting support from the motor. The so-called pedelec principle.

You can choose from five levels of electric support from the motor by means of the speed sensor alone, or in combination with the built-in throttle.

With a fully charged standard battery, you can reach distances of around 20-40 km, depending on how much you pedal, the support level you use, the topography and temperature, and it is possible to upgrade to a battery with a higher capacity or purchase an extra battery, if needed.

Gear system:

7 gears freewheel

7 gears with foot brakes and reverse possibility



Side-by-side P20



Side-by-side P20

Measures Side-by-Side:

Length	190 cm
Width:	100 cm
Boarding height:	17 cm
Wheel size:	20"
Max load:	210 kg

Color options (standard)

	Metallic red
	Metallic blue
	Metallic silver
	Matte black



m o b i l i t y

Bjerregårdvej 15, Timring · DK-7480 Vildbjerg · Telephone: +45 99 92 06 00 · www.pfmobility.dk