

EASYUP®
support arms

STAND-UP
SUPPORT ARMS



STAND-UP aid for at the toilet
Living independently for longer



Dutch quality

EASYUP®
support arms

The EasyUP STAND-UP support arms

A high-quality product that perfectly suits both home and healthcare situation, with the aim of increasing independence and preventing falls during every toilet visit.

Characterized by the combination of high-quality materials, easy installation, and lasting comfort, the EasyUP toilet support arms represent a reliable and user-friendly choice for every end-user.



Technical information


Constructed from
stainless steel

**Combinable with
wash-dry systems**

Support arms are
foldable and fixable

Supports up to 150 kg

Tested according
to ISO 17966



(A) Width between support arms 515 mm

(B) Suitable for mounting on toilets
up to 39 centimeters wide

(C) Support arm height adjustable in three positions
(Height from toilet seat: 25 cm, 29 cm, 32.5 cm)

(D) Can be attached to toilets with seat mounting
from 125 mm up to 200 mm (up to M8)

(E) Minimum distance from wall or vertical rim to
mounting holes up to 30 mm

EASYUP®
support arms

Onderdeel van Aerolet®

AEROLET®
Het Lentfert 38
7547 SR Enschede
The Netherlands
T: +31 (0)53 - 303 30 24
E: info@aerolet.nl

→
SCAN THE QR AND VISIT OUR
WEBSITE FOR MORE INFORMATION

