




























# Grillo A.Seating koot ja painot




	A	B	C1	C2	D	E	F	G	L	M	Käyttäjän pituus (H)	kokonais- paino	max. kuorm.
<b>mini</b>	16-26 cm 6,2-10,2 in	16-26 cm 6,2-10,2 in	30-41 cm 11,8-16,1 in	30-50 cm 11,8-19,6 in	12-15 cm 4,7-5,9 in	15 cm 5,9 in	49 cm 19,2 in	59 cm 23,2 in	15-25 cm 5,9-9,8 in	68 cm 26,7 in	75-105 cm 29,5-41,3 in	20,4kg 44,9 lb	25 kg 55,1 lb
<b>Small</b>	20-30 cm 7,8-11,8 in	20-30 cm 7,8-11,8 in	35-47 cm 13,7-18,5 in	35-60 cm 13,7-23,6 in	15-20 cm 5,9-7,8 in	17 cm 6,6 in	53 cm 20,8 in	65 cm 25,5 in	20-31 cm 7,8-12,2 in	68 cm 26,7 in	90-120 cm 35,4-47,2 in	22,6 kg 49,8 lb	35 kg 77,1 lb
<b>Medium</b>	25-35 cm 9,8-13,7 in	28-43 cm 11-16,9 in	40-55 cm 15,7-21,6 in	40-70 cm 15,7-27,5 in	16-20 cm 6,2-7,8 in	20 cm 7,8 in	58 cm 22,8 in	74 cm 29,1 in	25-40 cm 9,8-15,7 in	70 cm 27,5 in	110-145 cm 43,3-57 in	24,3 kg 53,5 lb	45 kg 99,2 lb

Viivoilla erotetut numerot osoittavat vähimmäis- ja enimmäissätöalueen.

Tessa ja hänen tarinansa odottavat sinua Ormesan virallisella verkkosivustolla yhdessä muiden fantastisen maailmanlaajuisen yhteisöimme päähenkilöiden kanssa.





**A special thanks to the families who participated with enthusiasm in our project. It is thanks to our customers and their stories that we find the energy and stimuli to evolve day by day.**

**CE** complies with MDR (EU) 2017/7

COMPANY WITH  
QUALITY SYSTEM  
CERTIFIED BY DNV GL  
= **ISO 13485** =

---

**Ormesa srl** - Via A. Da Sangallo, 1 | 06034 Foligno (PG) - ITALY  
Tel. +39 0742 22927 | FAX +39 0742 22637 | [info@ormesa.com](mailto:info@ormesa.com)

**[www.ormesa.com](http://www.ormesa.com)**