

ECO TRAVEL

THE SMART WAY TO TRAVEL



REHASENSE

**NEW
PRODUCT**



ECO TRAVEL

COMPACT & TRAVEL-FRIENDLY

Introducing the Rehasense ECO TRAVEL power add-on for manual wheelchairs—a game-changer for those seeking enhanced mobility on the go. Specifically designed with (air) travel in mind, this is the lightest power add-on in our lineup, offering unmatched convenience and portability. The ECO TRAVEL effortlessly transforms your manual wheelchair into a powerful, easy-to-maneuver companion, perfect for both short trips and extended journeys. Lightweight and portable – designed for easy air travel. It has high capacity battery, foot rest solution, weight, basket and tetra handles option. You can choose between the two versions ECO TRAVEL STANDARD or ECO TRAVEL ADVANCED.

EASY SETUP & PERFORMANCE:

- **Quick and Easy Installation:** The modular mounting system requires no pre-mounted brackets, allowing for hassle-free setup
- **Customizable comfort:** User-friendly settings, including height and length adjustments, ensure optimal sitting and driving comfort
- **Unmatched battery range:** Longest battery range on the market in that category – up to 45 km with std. battery
- **Enhanced Stability:** The 2-clamp attachment on the front frame provides optimal driving stability and protection

DOCKING & UNDOCKING

Swift Docking and Undocking: Experience the fastest setup on the market – quick and effortless, with only seconds needed to be ready to drive

SAFETY

- Dual braking system (drum and electronic brakes)
- Includes intelligent monitoring system to avoid any unfortunate situation
- The self-centering system makes driving safe and comfortable

COMPATIBILITY

- Can be used with rigid and folding manual wheelchairs - and also models with removable leg rests
- Can be mounted on almost all wheelchair models on the market*
- Covers mounting frame widths between 220-480 mm

COLOUR

Black



Green Orange Blue Red Black

*We do not recommend using ECO TRAVEL with carbon fibre wheelchairs. Rehasense takes no responsibility for any damage.

EXTRACT OF TECHNICAL DATA

	STANDARD	ADVANCED
Total weight	13 kg	14 kg
Transport weight	8 kg	9 kg
Max. user weight	110 kg	
Wheel size	8"/200 mm	
Engine power	250 W	350 W
Max. engine torque	19 Nm	24 Nm
Speed modes	2 gears	3 gears
Max speed	up to 15 km/h	up to 20 km/h
Climbing power	up to 10% (6°)	
Standard battery	Li-on, 48 V/ 5,8 Ah/ 278 Wh	
Driving range	up to 45 km	
Display	1.34" LCD black/white	1.81" LCD, color
• Setting options	None	Yes
• Language	English only	Multilingual
• Cruise control	No	Yes
• Motor brake (CLC)	No	Yes
• WALK function	No	Yes
Colors/Personalization	Black only	5 different
Handlebar	Fixed	Foldable
Handlebar stem	100 mm, fixed length	Length adjustable 180-215 mm and easily foldable (without tools)
Support legs type	Simple feet	Incl. wheels
Travel suitcase	Available as an accessory	Included

We reserve the right to make changes.



Rehasense's quality management system complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TUV NORD. Rehasense's medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR.) Rehasense's medical devices are UKCA-marked and meet the requirements of MHRA regulations.

