



YHTEENSOPIVA



Basmodellen för den förnyade Action-familjen

.....

Invacare **Action3 Basic** är ett av de bästa exemplen på Invacare-produkter anpassningsförmåga och kompatibilitet med Action-familjs tillbehör.

Justerbarheten och anpassningsförmågan hos **Action3 Basic** rullstolsdelar ger effektiviteten till produkt Distributörer och många alternativ för användarna.

Förnyade, fällbara, nya generationens **Action3 Basic** rullstolen kombinerar ny stilfull design, förbättrad användbarhet, stabilitet och komfort.



Yes, you can.®

Tekniska data


							
Drivhjul	405 / 430 / 455 / 480 mm	400 / 450 mm	0° / 3° / 6°	*(360 / 385 / 410) / 435 / 460 / 485 / 510 mm * mindre länkhjul	0° / 7° / 0° - 30° / -6° - +24°	350 - 510 mm	80° / 0 - 80° / 90°
Transithjul	405 / 430 / 455 / 480 mm	400 / 450 mm	0° / 3° / 6°	460 / 485 / 510 mm	0° / 7° / 0° - 30° / -6° - +24°	390 - 510 mm	80° / 0 - 80° / 90°
							
Drivhjul	200 - 400 / 330 - 460 / 440 - 510 mm	200 - 280 / 240 / 170 - 300 / 210 - 350 mm	Sitsbredd + 190 mm	295 mm	980 - 1105 mm	835 - 960 mm	820 - 1050 mm
Transithjul	200 - 400 / 330 - 460 / 440 - 510 mm	200 - 280 / 240 / 170 - 300 / 210 - 350 mm	Sitsbredd + 160 mm	260 mm	950 - 1050 mm	710 - 910 mm	920 - 1050 mm
							
Drivhjul	9°	870 mm	125 kg	14,2 kg	7,6 kg		
Transithjul	9°	870 mm	125 kg	13,5 kg	10,7 kg		

All product dimensions are in millimetres to apply EN 12183 standards. Our wheelchairs comply with ISO 7176-19 standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible. For more pre-sales information about this wheelchair, please consult the www.invacare.eu.com website. It can be found under the section Medium Active Wheelchairs / Action3NG / Sales and Images in the product catalogue. For information about typical user description, intended use, adjustments and transportation, please refer to the user guide.

Ramfärg



Hypnotic Purple



Matt Black




Blue Pop



Matt Red



Black Nylon
(klädsel)



Pearl White



Slate Grey



Electric Green



Anthracite Black

Tillverkare:

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196
Gland Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeuropa@invacare.com
www.invacare.eu.com

© Camp Mobility Finland 2019

