

Foot Control



Camp Mobility Foot Control for electric wheelchairs - New technical innovation!

Camp Mobility Foot Control for electric wheelchairs is widely patented technical innovation that works smoothly also in curves and tight turns. Easily manoeuvrable foot control works like a normal joystick - just with your ankle movements. Move ankle forward to go on, move backward to go back. Release center to brake, move side to turn around. - Easy!

Compatible with Universal R-Net, LiNX- and DX2-electronics.

Comes with additional button controllers that can be mounted to the back or where ever the user wants to. Foot Control can be easily changed to normal control or attendant control just pressing the button.

CAMP[®]
MOBILITY

Features



Works like a joystick - only with your ankle!

Move ankle forward to go on, move backward to go back. Release center to brake, move side to turn around.



ETC Patented - Compatible with R-net, LiNX- and DX2-electronics

Fits to the most electric wheelchairs with Universal R-Net, LiNX- or DX2-electronics.



Most of the basic features can be operated with Foot Control

Power ON / OFF, electrical adjustments and driving speed.



Attachment to the wheelchair's legrest

Stable and safe pedal works like a normal joystick, only with small ankle movements.



Comes with additional button controllers

that can be mounted to the back or where ever the user wants to. Foot Control can be easily changed to normal control or attendant control just pressing the button.

Technical data



Measurements:
260 x 165 x 50 mm



Weight: 1,25 kg

Sales Finland
Camp Mobility
Patamäenkatu 5
33900 Tampere
Finland

We preserve rights to make
changes to the product
specifications without
further notice
Camp Mobility © 2020

Want to become a dealer?
Contact:
info@campmobility.fi