

• MADE IN ITALY SINCE 1980 •

STANDY





A smart stander for your well-being.

Standy by Ormesa is a standing frame for adults and teens with various levels of ability and autonomy, and it is able to ensure an excellent verticalization.

Safe, versatile, easy to adjust: it provides **stability and comfort in every situation**, thanks to a set of accessories designed to meet the particular needs of each user.





Ergonomic

Innovative and functional design



The use of **Standy** facilitates the physiological verticalization of the spine, improving head control and horizontal gaze.











Technical features and adjustments



Frame

no edge frame, in order to enable trunk rotation and side-arm motion. Height adjustable frame with "safety locking system".



Tools free

adjustment without tools but with special knobs, designed by Ormesa, with effective and safe grip. The knobs effectively maintain tightening thanks to an internal mechanism that does not mark and ruin the tube.





Padded thoracic support

adjustable in height.





Tray


removable and adjustable in height and depth.





Padded back support

adjustable in depth with safety lock.



Padded knee pads

adjustable in width, depth and height.



Heel stops adjustable in width and depth.



Anti-slip footrest

with wheels for moving without user.



Components



865 Headrest adjustable in height, tilt, back and forwards.



857 Side support for the pelvis



945 Handle to open the back support



861 4 wheel set with brake



902 Front footrests with straps

Standy sizes and weight

A: length encumbrance

B: base length encumbrance

C: footrest to back support height

D: footrest to thoracic support height

E: base width encumbrance

F: footrest to knee pads height

G: width encumbrance

L: thoracic support width

H: user height

Total weight

Maximum load

standy 3

89 cm

76 cm

from 63 to 93 cm from 81 to 121 cm

69 cm

from 24 to 49 cm

74 cm 30 cm

from 110 to 160 cm

30 kg 90 kg standy 4

100 cm 85 cm

from 72 to 106 cm

from 84 to 129 cm


79 cm


from 30 to 57 cm

80 cm 39 cm

from 135 to 180 cm

38 kg 120 kg







COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 13485 =

Ormesa srl - Via A. Da Sangallo, 1 | 06034 Foligno (PG) - ITALY Tel. +39 0742 22927 | FAX +39 0742 22637 | info@ormesa.com

www.ormesa.com